

# **Olympus Spinning Disk - Widefield mode**

Written By: Joana Raquel Delgado Martins



# Step 1 — Widefield mode

Adjust Display		? 🕂 🗙	My Functions 1: Widefield	? 7
Histogram		▲ ▲	s Widefield	
🔍 10 bit	• 🕀 🖉 🔛 😣			
A			SF open	® Transmitted
	the star star star star		🔰 🔍 dapi wf	V DAPI/Cy3
< 300 400 S	500 600 700 800 900 10	>		
Min:		Max:	GFP WF	Reflected 2
Mean Intensity:	258.64	2825		
Pixel Count:	4'186'574		Cy3 WF	
OFixed Scaling				
Left:	Right:			
104	1183		Cy5 WF	
Auto Contrast				
Left:	Right:		c	>
0.1 %	0.1 %		•	

- To use the microscope on the Widefield mode different observation methods are available:
- You can find them under My functions - "Widefield" for easier access during Live mode.
  - DAPI WF
  - GFP WF
  - CY3 WF
  - Cy5 WF
  - DAPI/Cy3 (using the dual camera mode)
  - Be aware that the Cy5 WF observation method might also show signals is the red range as a broad light source is used (no excitation filter is present).

# Step 2 — Design Experiment

W4         W4         W         Experiment Manager         7 • ×           Recogan dia tones         Image: Comparison of tones         Image: Comparison of tones         Image: Comparison of tones		🖣 💽 Start Page 🗙 💽 Experiment* 🗙 🕒 dapi 488 561 test timelaps	e.oex*
Excide system         Image: Control of the system           Appr         Dolluit           New Experiment         srt           Stop -	Eventer Templates     Cont ZDC with MultiPos Sem     System Templates     Cont ZDC with MultiPos Sem     Cont ZDC with MultiPos Sem	¥•⊛≋≍⊘ ⋬∑≧∎∙ ⊛• የ•¢•≪•●• ∦•	
Display Enhancement Experiment Bightone: 10 12 Experiment (2)	Add Current Experiment to Favorite Templates Add Current Selection to Favorite Templates Widefect DAPI CFF Cy3 Cy5		
Contrast:         D0 (c0.00)           50 (c)         Obligation:           University         Obligation:			
Append the lapse loops     Count      C			
▲ LeartID Controler ■ 645 cm ■ 0 0 0 5 %		DAPI WF GFP WF Cy3 WF Cy5 WF	
		20.02 ms 20.02 ms 20.02 ms 20.02 ms	
91m 1 9 3 %			
Speriment Manager			

- Start a new Experiment on the "Experiment manager" tab.
- For a 4 color sequential image you can use the template Favorite Templates "Widefield DAPI GFP Cy3 Cy5".
  - You can also design your experiment using the available widefield observation methods from the menu.
  - For multichannel images do not forget to add the respective option.

### Step 3 — Adjust LED intensity and exposure time

Camera Control 7 * X AF • • 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2	Histoscope Control 7 9 X Histoscope Control 7 9 X Exclude spikes in histogram
Lve Snapshot Movie recording Snapshot Movie recording Subtrance Live Snapshot Sn	Apply Default Display Enhancement Sol 50 Contrast: 50 Con
1  Resolution: 2048 x 2048  Live/Novie: 2049 x 2048  Snapshol/Process: 2049 x 2048  V	Carrents: 1 2 2 2 2 2 2 2 2 2 2 2 2 2
Apply to Live	20 µm ▲ Lasert.ED Combiner 405 nm 43 \$% 445 nm 0 \$% 445 nm 40 \$% 458 nm 40 \$% 103 Aceture 100 ↓ 103 Aceture 100 ↓ 103 Aceture
Stage Navigator 7.9 ×	S ia mm 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0

- Double click on the observation method box to activate it.
- Check the "Synchronize Shutter" box option if you would like to minimize exposure of your sample to light.
- Press Live.
- Adjust exposure time accordingly.
- You can also change LED intensity under the "Microscope Control tab":
  - Scroll down until you find pE-300

(i) the appropriate LED will be active and can be changed for each channel.

# Step 4 — Check image histogram and avoid clipping



A Make sure your image is not saturated with the help of the Adjust Display pannel.

- Here you can toggle the camera bit-depth between 16 and 8-bit .
- Using the camera at 16-bit (recommended), Max intensity should be less than 65 536. At 8-bit, the Max intensity should stay bellow 255.
- For more details refer to the different camera modes here: <u>Olympus Spinning Disk 2:</u> <u>Multichannel + Z-stack acquisition</u>.

# Step 5 — Get settings



- Once you are satisfied with your settings right click on the observation method and select "Get settings".
- Repeat the procedure for the remaining observation methods.
- Name your experiment.
- Press Start.
- Visualize your image