


# Leica SP8 upright BL3 (Irchel) - 1: Start-up

How to start up and mounting sample on Leica SP8 upright confocal laser scanning microscope located at Irchel Campus, at the Medical Virology (Y36-M-92).

Written By: Caroline Aemisegger




# INTRODUCTION

How to start up and mount your sample on the Leica SP8 upright confocal laser scanning microscope located at Irchel Campus in the Medical Virology (Y36-M-92).


Please find detailed information about the system setup [here](#).

## Step 1 — Switching ON hardware




- Switch ON the **fluorescence lamp**.
  - Switch ON the "**Scanner Power**", "**Laser Power**", and turn the "**Laser Emission**" key to "On-1" (control unit underneath the table).
  - Switch ON the **power knob** (on the PC table).
- ❗ *The PC and microscope will be switched on.*

## Step 2 — Sign-in




- **Sign-in** with your **ZMB core** credentials.

## Step 3 — Starting up the "LAS X" Software



- **Start** the "LAS X" software.
- Make sure "**machine.xlhw**" is selected as "Configuration", and "**DM6000**" as "Microscope".
- Select either "Resonant" (ON) or non-"Resonant" (OFF) mode.
- Click "**OK**".
- Click "Yes" in order to **initialize the x/y stage**. *Please make sure nothing is placed currently on the stage.*


## Step 4 — Switch ON the lasers in the software



- Go to **"Configuration"**.
- Select **"Laser Config"**.
- **Switch ON** the lasers you will need.
- Go back to **"Acquire"**.





## Step 5 — Choosing an objective




- Lower the stage by pressing the "LOWER Z" button on the right side of the microscope.
  - ⚠ This is a mandatory step as it avoids possible collision of the stage and objective during exchange of inserts and/or samples.
- You can now toggle between objectives within the software (drop-down menu).
- Select the **10x dry objective** .
  - ⓘ In order to facilitate the focusing it is recommended to start with the 10x dry objective.

## Step 6 — Mount and position your sample



- **Insert** your sample with the **coverslip facing up** and fix it with the two springs.
- Move your sample under the objective with the help of the external controller "**Smart Move**".
  - Movement in y-direction.
  - Movement in x-direction.
  - Toggle between coarse movement "XY Fast" and slow movement "XY Precise".

## Step 7 — Choose the Eyepiece light path



- On the touch screen at the microscope stand choose the **objective tab**.
- Make sure the **eyepiece** is selected as the "Port - Top Port 100%".




## Step 8 — Focus your sample




- On the touch screen at the microscope stand choose the **light path tab**.
- Click "FLUO" and choose an appropriate "FLUO-Filtercube" : e.g. "GFP".
- Open the "IL -Shutter" (if activated the dot is yellow).
- Look through the oculars and focus your sample by using:
  - the focus wheel on the microscope stand,
  - or the z-wheel on the external controller (Smart Move).
- ① Moving sample upwards (towards objectives) turn z-wheels clockwise/away from you. Moving sample downwards (away from objectives) turn z-wheels counter-clockwise/towards you.
- Toggle between "Z FINE" and "Z COARSE" directly on the Smart Move.

## Step 9 — Optional - Save your focus position




- i** The storage of the focal plane is helpful in order to find the focus back if the sample or objective will be changed.
- To save your current focus position select the "**xyz tab**" and the "**Focusdrive Z**" on the touchscreen of the microscope.
  - Click the "**Upper Focus Limit**" button.
  - Press "**Set**".
    - If done successfully you will see an upper marker line appearing.
  - Press the "**Lower Limit**" button in order to move down (*for safe change of the objective or the sample*).

## Step 10 — Switching to a higher magnification



- Remove your sample and toggle within the software to the objective of choice.
- ❗ Depending on the objective different **immersion media** will be used. Apply directly on the sample.
  - Oil objectives: "Type-F" immersion liquid.
  - "Glycerin" objectives: "Type-G" immersion liquid.
  - "Water" objectives: ddH<sub>2</sub>O (always use fresh).
- ⚠ Please further consider the additional information in the next step to guaranty proper image acquisition.
- **Mount** your sample again and press the "**Upper Focus Limit**" button.
- Focus your sample as described previously.

## Step 11 — Additional information - Immersion objectives



⚠ For some objectives the **correction collar** has to be adjusted.

- **20x IMM** (multi-immersion - Oil, Glycerin or Water) needs to be set to the corresponding **immersion media** ("OIL", "GLYC" or "0.17-W" (with cover glass) or "W-0" (without cover glass)).
- Make sure that the cap of the **spring-loaded front lens** is released (working position).

⚠ Please, **DO NOT** remove the objectives for adjustment. They can be also accessed on the system.