

Leica SP5 Mid UV-Vis (Irchel) - Part 1: Start-up

How to start up and mounting sample at the SP5 Mid UV-VIS confocal laser scanning microscope located at Irchel, room Y34-E-36.

Written By: Jana Döhner



INTRODUCTION

How to start up and mount your first sample on the SP5 Mid UV-Vis confocal laser scanning microscope located at Irchel, room Y34-E-36.


Please find detailed information about the system setup [here](#).

Step 1 — Switching ON hardware



- Switch ON the **fluorescence lamp**.
- Switch ON the "**PC /Microscope**", "**Scanner Power**" and "**Laser Power**" and turn the "**Laser Emission**" key to "On-1" (main switch board underneath table right hand side).

Step 2 — Sign-In




- **Sign-in** with your **ZMB core** credentials.

Step 3 — Starting up the "LAS AF" software



- Start the **"LAS AF"** software.
 - If the **"Resonant Scanner"** is needed check "Activate Resonant Scanner".
 - Click **"OK"**.
 - Select **"Yes"** in order to **initialize the x/y stage**. *Please make sure that nothing is placed on the stage.*
- ① An x/y stage initialization is necessary in order to use the "Tilscanning", and "Mark and Find" function.

Step 4 — Switch ON the lasers in the software



- Go to "**Configuration**".
- Select "**Laser**".
- **Switch ON** the lasers you will need.
 - Adjust the **Argon** laser to **20%**.
 - The **405 nm** laser has to be switched ON via the **external** button.
- ① Follow next step, if in addition the **355 nm laser** is needed, otherwise continue with step 6.


Step 5 — Optional - 355 nm UV-Laser



i **Optional** step: follow this step only if the 355 nm laser is needed for your experiment, otherwise continue with step 6.


- Start the 355 nm laser.
 - Turn ON **main switch** (on the back). The LED "SYSTEM FAULT" lights up. *Please note, key switch in the front must be off.*
 - Wait until laser has powered up and the indicator "SYSTEM FAULT" has turned off. *This takes approx. 1min.*
 - The "INTLK OK" should light up now.
 - Turn the **laser key** to "ON".
 - **Set** the power with the control knob. *Turn for fine adjustment, press and turn for coarse adjustment. The LED "LASER ON" lights up.*

Step 6 — Choosing an objective




- Go back to **"Acquire"**.
- **Lower the objective turret** by pressing the downwards "Z" button on the right side of the microscope.
 - ⚠ This is a **mandatory step** as it avoids possible collision of the objectives and stage during exchange of inserts and/or samples.
- You can now toggle between objectives within the software (drop-down menu).
- Select the **10x dry objective**.
 - 📌 *In order to facilitate the focusing process it is recommended to start with the 10x dry objective.*

Step 7 — Mount and position your sample




- **Push** the condenser arm of the microscope to the **back**.
- Insert your sample with the **coverslip facing down** and fix it with the two springs.
 - ❗ *Special stage inserts/adapters are available for other samples than regular slides (please see last step of this guide).*
- Move your sample over the objective with the help of the external controller **"Smart Move"**.
 - Movement in y-direction.
 - Movement in x-direction.
 - Toggle between coarse movement "XY Fast" and slow movement "XY Precise".
- ☑ Bring back condenser arm to its straight position.

Step 8 — Focus your sample




- On the front panel of the microscope select an appropriate **fluorescence filter**:
 - filter for UV dyes like DAPI,
 - filter for green dyes like FITC or Alexa 488,
 - filter for red dyes like TRITC or Alexa 568.
- Press the **"SHUTTER"** in order to illuminate your sample.
- Look through the oculars and **focus** your **sample** by using the **focus wheel** on the microscope stand or the **external controller (Smart Move)**.
 - ① Moving objectives upwards (towards sample) turn z-wheels clockwise/away from you. Moving objectives downwards (away from sample) turn z-wheels counter-clockwise/towards you.
- Toggle between "Z FINE" and "Z COARSE" directly on the Smart Move.

Step 9 — Optional - Save your focus position



- ❗ The storage of the focal plane is helpful in order to find the focus back if the sample or objective will be changed.
- **Jointly** press the "SET" and upper "Z" button (right side of the microscope stand) in order to set the current z-positon to zero.
- ❗ *Depending if already a focus was saved by a previous user, you have to do that "once" (nothing was saved) or "twice" (different focus was saved and needs to be first deleted).*
- The z-position on the display should now show "**0 mm**".
- Press the lower "**Z**" button in order to move down *(for safe change of the objective or the sample)*.

Step 10 — Switching to a higher magnification



- Remove your sample and toggle within the software to the objective of choice.
- ❗ Depending on the objective different **immersion media** will be used. Apply either on the sample or directly to the objective.
 - Oil objectives: "Type-F" immersion liquid.
 - "Glycerin" objectives: either "Type-G" immersion liquid (for RT measurements) or "Glycerin" immersion liquid (for measurements at 37°C).
 - "Water" objectives: ddH₂O (always use fresh).
- ⚠ Please further consider the additional information in the next step to guaranty proper image acquisition.
- **Mount** your sample again and press the upper "Z" button.
- Focus your sample as described previously.

Step 11 — Additional information - Immersion objectives





⚠ For some objectives the **correction collar** has to be adjusted.

- **20x IMM** (multi-immersion - Oil, Glycerin or Water) needs to be set to the corresponding **immersion media** ("OIL", "GLYC" or "0.17-W" (with cover glass) or "W-0" (without cover glass)).
- **40x OIL** fully open the aperture (NA 1.25) of the lens by turning the correction collar in clockwise direction.
- You can correct for the **cover glass thickness** and temperature at the **63x Glyc** (0.14-0.20) and **63x water** (0.14-0.18). Standard is usually 0.17 mm
- Make sure that the cap of the **spring-loaded front lens** is released (working position).

⚠ Please, **DO NOT** remove the objectives for adjustment. They can be also accessed on the system.

Step 12 — Additional information - Stage inserts



i Stage inserts can be easily exchanged or adapted for other samples than regular slides. You can find additional inserts and adapters in the box on the table behind the microscope.

- Standard insert for slides, including chamber slides.
- Standard insert with petri dish adapter (for 36 mm diameter).
- Insert for multiwell plates.