

Leica SP8 Falcon (Irchel) - Part 1: Start-up

How to start up and mount your first sample on the Leica SP8 FALCON confocal laser scanning microscope located at the Irchel Campus, room Y42-H-81.

Written By: Joana Raquel Delgado Martins



INTRODUCTION

How to start up and mount your first sample on the Leica SP8 Falcon confocal laser scanning microscope located at the Irchel Campus, room Y42-H-81.

Please find detailed information about the system setup here.

Step 1 — Switching ON Hardware



Switch on the red button (underneath the table on the left).

(i) Turns on fluorescence lamp and monitors.

- On the right side of the table:
 - Switch ON the "PC/Microscope", "Scanner Power" and "Laser Power" switches.
 - Turn the "Laser Emission" key to "ON-1".

Step 2 — Sign-In



• Sign-in with your ZMB core credentials.

Step 3 — Starting up the "LAS X" Software

Feedback Citric Receiver LAS X	Leica Application Suite X 3.5.6.21594 Configuration : Microscope : DMB =	Leica Application Suite X 3.5.19976 Microscope Stand
Matriacree 1 for LAS X Matriacree 2 for LAS X	Resonant : Load settings at startup : Copyright © 2019 Leica Microsystems CMS GinbH	Initialize Stage: DMIB Stage? If yes, please protect microscope condenser by limiting stage working area! Bi Dialog will close automaticaly in 00:01:00 with No selected.

- Start the "LAS X" software. Select:
 - "machine.xlhw" for "Configuration",
 - and "DMI8" as "Microscope".
- Select either "Resonant" (ON) or non-"Resonant" (OFF) scanning mode.
 - *Use "Resonant" scanner for fast acquisition and/ or live imaging. However, not advised for FLIM!*
- Click "OK".
- Click "Yes" to initialize the x/y stage. Make sure nothing is placed on the stage.
 - (i) An x/y stage initialization is necessary to use the Navigator function.

Step 4 — Switch ON the lasers in the software

Retree	TCS SP8 🗘 Configuration Acquire	
Beam Path		Currently available Lasers Adjust Laser Settings *
USB Panel Objective		Diode 405 : Standby Image: 440pulsed : Standby Image: Will : ON Image: Standby
Hardware		

- Go to "Configuration".
- Select "Laser Config".
- Switch "ON" the lasers you will need.
 - When "ON", the WLL should be at 85% by default.
- Go back to "Acquire".

Step 5 — Choosing an objective



 Lower the objective turret by pressing the "Z downwards" button on the right side of the microscope.

A This step avoids possible collision during placing of inserts and/or samples.

- Select the **10x dry objective** via the "LAS X" software.
- (i) In order to facilitate the focusing it is recommended to start with the 10x dry objective.

Step 6 — Check stage insert



- Choose the appropriate sample holder:
 - The depicted stage insert is usually placed at the microscope.
 - Other holders (eg. 96 well plates) are stored in a box on a shelve behind the microscope.

(*i*) If necessary, stage inserts can be easily exchanged as they are held by magnets.

Step 7 — Mount and position your sample



- **Push** the condensor arm of the microscope to the **back**.
- Insert your sample with the **coverslip facing down** and fix it with the two springs.
- Move your sample above the objective with the help of the external controller "Smart Move".
 - Movement in <u>y-direction</u>.
 - Movement in <u>x-direction</u>.
 - Toggle between coarse movement "XY Fast" and slow movement "XY Precise".

Bring back condenser arm to its straight position.

Step 8 — Optional - Exchanging stage inserts



- To exchange a stage insert **Push** the condensor arm of the microscope to the **back**.
- Remove the stage insert and place the appropriate one.

Step 9 — Focus your sample



- On the touch screen at the microscope stand choose the light path tab.
- Click "FLUO" and choose an appropriate "FLUO-Filtercube" : e.g. "GFP".
- Open the "IL -Shutter" (if activated the dot is yellow).
- Look through the oculars and focus your sample by using:
 - the focus wheel on the microscope stand,
 - or the z-wheel on the external controller ("Smart Move").
 - (i) <u>Turn z-wheels clockwise</u> to move objectives upwards (closer to the sample). <u>Turn z-wheels</u> <u>counter-clockwise to</u> move objectives downwards (away from sample).
 - Toggle between "Z FINE" and "Z COARSE" directly on the Smart Move.

Step 10 — Switching to a higher magnification



- Remove your sample and toggle within the software to the objective of choice.
- Depending on the objective different immersion media will be used.
 Apply directly on the sample.
 - <u>Oil objectives</u>: "Type-F" immersion liquid.
 - <u>Glycerin objectives</u>: "Type-G" immersion liquid . For room temperature use the 23°C glycerin media. For live cell imaging the 37°C one.
 - <u>Water objectives</u> : Use fresh double destiled water.
- You can move (back and forth) the condenser arm for ease of access.
- Please consider the additional information in the next step to guaranty proper image acquisition.
- Focus your sample as described previously.

Step 11 — Optional - Save your focus position



- (i) The storage of the focal plane is helpful in order to find the focus back if the sample or objective will be changed.
 - To save your current focus position select the "xyz tab" and the "Focusdrive Z" on the touchscreen of the microscope.
 - Click the "Upper Focus Limit" button.
 - Press "Set".
 - If done successfully you will see an upper marker line appearing.
 - Press the "Lower Limit" button in order to move down (for safe change of the objective or the sample).

Step 12 — Additional information - Immersion objectives



A For some objectives the **correction collar** has to be adjusted.

- 20x IMM (multi-immersion Oil, Glycerin or Water) needs to be set to the corresponding immersion media ("OIL", "GLYC" or "0.17-W" (with cover glass) or "W-0" (without cover glass)).
- (i) 40x water and 63x glycerol you can correct for the cover glass thickness (0.14-0.18 mm). Standard is usually 0.17 mm.
 - 40x water you can adjust for the correct cover glass thickness.
 - 63x glycerol you can adjust for the cover glass thickness of the corresponding temperature.
 Upper row for 23°C with the indicated 0.17mm.
 - Lower row for 37°C and indicated 0.17mm.
- Make sure that the cap of the spring-loaded front lens is released (working position). <u>Mandatory</u> for <u>all</u> immersion objectives.

Nease, DO NOT remove the objectives for adjustment. They can be also accessed on the system.