


TEM - FEI Tecnai Spirit

How to start up and image at the TEM - Tecnai Spirit located at the Center for Microscopy and Image Analysis (ZMB), Irchel Campus, room Y42-H-95, Zurich.

Written By: z mbstaff




INTRODUCTION

How to start up and image using the TEM - Tecnai Spirit located at the Center for Microscopy and Image Analysis (ZMB), Irchel Campus, room Y42-H-95, UZH, Zurich, Switzerland.

Our Tecnai Spirit TEM is an easy to use 120 kV transmission electron microscope (TEM) designed to provide high-contrast, high-resolution imaging and analysis in life science applications.

Step 1 — TEM - FEI Tecnai Spirit



- Apertures
- Goniometer (Stage)
- Liquid nitrogen tank
- Control software interface
- Phosphor screen
- Control panels
- EDX-detector

Step 2 — Refill the liquid nitrogen



- Fill the liquid nitrogen dewar flask from the tank in front of the facility.

⚠ Wear the supplied safety gloves and glasses, when you work with liquid nitrogen.

Step 3 — Fill the cooling system




- Wear safety glasses.
- ❗ Check the nitrogen level during your session. One filling lasts for up to two hours.

Step 4




- ❗ Turn on screens if not already on.
- The control computer of the Spirit is never logged of.
- ❗ In case you need to log in: User: **Supervisor** / PW: **supervisor**
- ⚠ Please, never log off Spirit computer.

Step 5 — Software on the microscope




- ★ Make sure the following softwares are running:
 - TEM User Interface
 - TEM Imaging and Analysis
 - Gatan Digital Micrograph
- Remote desktop connections can be set up with these icons:
 - Support PC
 - EDX PC

Step 6 — Login at Support PC



- Click on the link to open a remote desktop connection to the support PC
- ❗ For security reasons, the FEI Tecnai Spirit is not directly connected to the network, therefore you need to first save your images to the **"SpiritData"** folder. This drive is shared with the Support PC and you must transfer your images to your data drive on **"Z"** at the end of your session.
- Open the **"SpiritData"** folder in your account. And in a second window open your ZMB folder: **Z:Data/Your Folder**
- In the **"SpiritData"** drive create a folder with your name. This is the location where you have to save your data.

Step 7 — Checking the vacuum condition.



- Check if the vacuum is highlighted in green.
- ① "Gun/Col" should be around Log 6
- ① "Camera" should be around Log 26
- Turn on "High Tension". (Yellow = On)
- Turn on "Filament" (Yellow = On)

Step 8 — Aperture settings for high resolution brightfield imaging




- The condensor aperture needs to be inserted.
- The objective aperture needs to be inserted .
- The SA aperture needs to be retracted for TEM bright field imaging.
- The condensor and the objective aperture are both in position three for regular imaging .
- ① The position is indicated by a small pin.

Step 9 — Adjustments for alignment



- Set the microscope to Spot size 1 and to "Microprobe"
- This list shows the functions of the keys on the control panels.
- Magnification.
- Adjust the magnification to "**2550x**" for the alignment.
- ① Right control panel.
- Press the button L3 to normalize all lenses.
- ① Left control panel.

Step 10 — Beamshift alignment




- Use the trackball for quick beam shift adjustment.
- Intensity controller.
- ① "Fine" and "Course" change the sensitivity of the knob
- Click on "Col. Valves Closed" to open the column valves.

Step 11 — Beamshift alignment



- Turn the intensity controller either left, or right until you see the electron beam on the phosphor screen.
- Use the intensity controller to change the diameter of the beam. Make it as small as possible.
- Use the trackball to center the focussed beam in the middle of the screen.

Step 12 — Beamshift alignment



- Increase the magnification and check if the beam is still centered. If not center it with the trackball
 - Magnification.
- ❗ Check the centering of the beam up to 16500x
 - Trackball. (left control panel)
 - Go back to 2550x magnification.
- Turn the intensity button (left control panel) clockwise until the the screen current in the software is 2.3 nA or lower.
 - Screen current.

Step 13 — Check the condensor aperture alignment




- The beam needs to stay centered, if you defocus it by turning the intensity knob clock wise.
- If the beam does not stay centered, you have to align the condensor aperture.

Step 14 — Align the condensor aperture





- Center the condensor aperture by turning its X and Y axis adjustment.
- Center the beam.
- Centered condensor aperture.

Step 15 — Prepare the stage for sample insertion



- Activate "Int. Zoom", if it is not activated.
- Close the column valves.
- Go to the menu "**Stage2**".
- Expand the window.
- Go to the tab "Control".
- Click reset the "Holder"

Step 16 — Insert sample into TEM holder.




① Use non magnetic forceps. Use a magnification glass for better visibility if needed.

⚠ Do not touch the front part of the holder with bare hands.

① Everything in front of the sealing O-ring goes into the vacuum.



- Open the holder clamp with the supplied needle.
- Place the grid into the groove with the forceps.
- You can position the grid by gentle tapping the finger against the transparent tube of the TEM holder.
- Close holder clamp with the needle!

Step 17 — Insert the holder




- Align the small pin to the arrow on the stage.
- Insert the holder as straight as possible, until you feel a resistance.

Step 18 — Insert the holder




- If this message appears: Select the Single Tilt Holder
 - And confirm.
- Wait for the countdown to finish.
- Wait until the red light turns off.

Step 19 — Insert the holder



- Turn the holder counter clockwise, until it is dragged in to the machine.
- Make sure the holder is completely inserted.
- ① Guide the holder with your hand during insertion.

Step 20 — Start the live view




- Open the column valves.
- Go to 2550x magnification.
- Check if the screen current is still at 2.3 nA or lower.

⚠ Do not insert the camera if the screen current is higher than **2.3 nA, as this can damage the camera chip.**

- Go to "Digital Micrograph" on the right screen.
- Set the exposure time for the live view to 0.1 seconds.
- Click on Start view.
- Click on "Yes" to insert the camera.

Step 21 — Adjust the illumination




- Always check, that the intensity is in the green or yellow range.

⚠ To much intensity can damage the camera chip!

- If you see this "honeycomb" like structure, the intensity is by far to high and you need immediately to dim the electron beam, or retract the camera out of the beam.
- Digital Micrograph can show you an intensity value for each point in an image.
 - Move the mouse to a white area in your live view.
- ① The value for a white area should be around 7000 or lower. Adjust the intensity if that is not the case.

Step 22 — Adjusting the Z-height



i The microscope performs best, if the sample is set close to eucentric height.


- Press Eucentric focus on the right panel to reset the focus of the lenses.
- Press L2 (Alpha Wobbler) on the left control panel. The stage begins to tilt between +/-15°
- Reduce the movement of the image in the live view of Digital Micrograph.
 - by changing the Z height with the Z axis button on the control panel
- sensitivity can be adjusted by + and - buttons above the joystick
- The Z-height is shown in the User Interface on the left screen. Use it as an orientation. (Z-values are usually between -20 μm and +20 μm)

Step 23 — Video Alpha wobbler




- Reduce the movement of the image in the live view of Digital Micrograph.

Step 24 — Fine focussing with electronic wobbler




- Push the "Wobbler" button on the right control panel.
- Turn on the "Focus" knob, until you do not see a double image anymore.

Step 25 — Activate Live FFT (optional)




- i The Live FFT can be used to judge the focus and astigmatism of an image.
- Choose the "rectangular ROI" tool.
- Draw a square shaped region of interest with the mouse. You need to press the alt key, while drawing the region.
- Go to "Process", choose "Live" and click on "FFT".
- The live Fourier Transformation will appear in a separate window.
- i The image is in optimal focus, when the white disk is as wide as possible and no additional rings are visible.

Step 26 — Correct objective lens astigmatism with Live FFT




- To check for astigmatism go out of focus (overfocus), until you see rings in the FFT.
- If the rings are oval and not round, you need to correct for astigmatism.
- Choose "Stigmator" in the menu in TUI.
- Click on "Objective".
- Use the multifunction knobs on the left and right control boards to adjust the astigmatism. Check in the Live-FFT, if the rings get circular.
- Click on "None", when you are finished with the adjustment.
- ① It is possible to choose between three different presets. You can also duplicate presets by a right mouse click to have a setting as a backup. You can then get back to the original astigmatism settings, if your corrections made it worse.

Step 27 — Acquire an image with Gatan Digital Micrograph




- The intensity value should be in the lower green range.
- The setup should be set to "Record". Choose the "Exposure Mode" "User" and an exposure time of 0.5 sec.
- Click on "Start Acquire" to take the image.

Step 28 — Save an image in Gatan Digital Micrograph



- Save your image as TIFF.
- Or save your image as DM3.
- ① DM3 allows for more adjustment in the post processing, but you will need the Gatan software to open it.


Step 29 — Set up autosaving



i There is a way to automatize the naming and saving of acquired images.


- Go to the "ZMB" tab.
- Choose "Prepare for save as numbered".
- Follow the setup by typing in the storage location, name and starting number.
- After you prepared the saving path and name you can choose one of the auto numbered options (option 2, 3 or 5) to save new images.

Step 30 — How to use tracking and save points of interest.




- Open the "Stage" tab.
- Activate "Tracks" to visualize where on your sample you already have been.
- Click on "Add" to save points of interest.
- If you click on "Go" you can move the stage back to the saved coordinates.
- With a click on the magnifying glass you can see the tracking in more detail.

Step 31 — How to transfer my data?




- Go to the "SpiritData" folder on the microscope control computer.
- Create a new folder with your name.
 - Save all your images in this folder.
- Open the remote desktop connection to the support PC.
 - Copy your files from the "SpiritData"
 - to your data drive in Z.

Step 32 — How to end my session?



- Log off from the Support PC, when you are finished with your data transfer

Step 33 — Take out the holder



- Reset the stage.
- Close the column valves!
 - In this position the holder can stay, if you need to reposition your hands.
- Pull the holder out completely.

⚠ Do not tilt the holder during removal or insertion. An angle can lead to a vacuum failure of the system.

The Tecnai™ Spirit TEM is an easy to use 20 kV to 120 kV transmission electron microscope (TEM) designed to provide high-contrast, high-resolution imaging and analysis in life science applications.