

Olympus Spinning Disk - DIC images

Written By: Joana Raquel Delgado Martins



Step 1 — DIC observation method

N-4 C/P Experiment Manager 7.9 X Anappen of all times Comment Manager 7.9 X		M1-43 🗳 FM1-43 💐 488_RAPP 💐 RAPP_OUT 💐 405/640 🧩 🦉	?QQQ\$} ?
Exclusive protein in integrate Image: Constraint of the constr		Þ ♂ マ Adjust Display ? ♥ X 0 400 420 µm Histogram ▲ ^	My Functions 1: ZMB DIC ? # ×
Display Exhancement Control Control Control	 Current Settings Current Settings ABN A	0 480 430 gr n notogram notogram 100 gr notog	

- Start a new Experiment on the "Experiment manager" tab.
- Make sure the microscope is properly aligned and that you have checked your Koehler illumination.



- Add your observation method DIC to the experiment.
- Do not forget to add an additional DIC_OUT observation method at the end to remove the extra components.
- An extra DIC_OUT image will be generated, you can delete or ignore it. For more info please contact the ZMB staff.
- You can find also find DIC and DIC_OUT observation methods under My functions "ZMB DIC".

Step 2 — Adjust exposure time and polarizer



- Double click on the observation method box to activate it.
- Check the "Synchronize Shutter" box option if you would like to minimize exposure of your sample to light.
- Press Live.
- Adjust exposure time accordingly.
- Under the "Microscope Control tab" your will find the two DIC Polarizers:
 - Scroll down until you find IX3 DIC slider.

(i) Here you can adjust the orientation of the polarizer.

Step 3 — Check image histogram and avoid clipping



A Make sure your image is not saturated with the help of the Adjust Display pannel.

(i) The Camera can be run at 8 bit, 12 bit as well as 16 bit (recommended).

- Here you can toggle the camera bit-depth between 16, and 8-bit and adjust the pixel clock (camera read out speed).
- ▲ If you are acquiring at 8 bit the histogram max should be lower than 255, at 12 bit lower than 4095, while if you are using 16 bit the max should be lower than 65535.
- For more details refer to the different camera modes here: <u>Olympus Spinning Disk 2:</u> <u>Multichannel + Z-stack acquisition</u>.

Step 4 — Get settings and start



- Once you are satisfied with your settings right click on the observation method and select "Get settings".
- Repeat the procedure for the remaining observation methods.
- Name your experiment.
- Press Start.
- Visualize your image

Step 5 — Fluorescence + DIC



If you want to acquire fluorescence + DIC images, make sure you remove DIC components from the light path after the DIC image and before acquiring the next fluorescent image, time point, position etc.

- Remove DIC components by adding a DIC_OUT observation method.
- (i) An extra DIC out image will be generated, you can delete or ignore it. For more info please contact the ZMB staff.
- Use the multicolor option to overlay the two channels.
- The 2 channels should have the same bit depth, size and format.