


Leica LAS X Navigator - Intro

This guide of the Center for Microscopy and Image Analysis shows the use of the Navigator on a Leica CLSM SP8 and Stellaris 5.


Written By: Karin Seubert



INTRODUCTION


This guide of the Center for Microscopy and Image Analysis shows the use of the Navigator on a Leica CLSM SP8 and Stellaris5 for getting an overview over the sample and setting different points for data acquisition. For the acquisition of tile scans, please refer to the guide ["Leica LAS X Navigator - Merge regular and irregular shapes"](#).

Step 1 — Starting the Navigator




- ❗ Set up your basic image acquisition for an overview. If you need a big overview (different tissue slices, well plates, etc.), consider using a low magnification dry objective (eg 10x) before applying immersion media.
- ❗ For the fastest overview, zoom out using the smart wheel on the panel bellow the computer monitor.
- ⚠ Make sure you initialized the stage during the startup. For more info refer to the the microscope startup guide.
- Start the **Navigator** by clicking on the symbol.
- A new window appears with the pannels accessible on the left.

Step 2 — Choose a Sample Carrier




- Go to the **"Carrier"** tab.
 - Choose your sample type from one of the drop-down menus.
- A preview of your sample holder with your sample (here, a coverslip on a slide) is shown.
- In the main window, you see a bigger version of it with a little square marking the position of your field of view.
- Gently move in x/y with the external controller **"Smart Move"** in **"XY Precise"** mode and observe the square moving.
- ⚠ Please be aware that this unaligned mode only gives you an idea where you are. Make sure to always stay in the indicated area to avoid crashing of the objectives.
- If you want to align your sample, click on **"Align"** and follow the instructions (best done with a 10x objective).

Step 3 — Getting an overview



- Use **"Fast Live"** to get a preview of the region you are currently imaging.
- Use the mouse left button and drag to navigate to a different stage position. With the mouse wheel you can zoom in and out.
- Double click (left mouse button) on the current position to get a snapshot of this position.
- Click on **"Spiral"** go get a quick overview of the sample with the field of view spiraling out.
- Press **"Stop"** once you are happy with the overview area.
- Choose which **channels** should be displayed using the check boxes.
 - Here you can **Auto-scale** the channels.
- Center your current field of view if you ever get lost.


Step 4 — Define your areas of interest



- Select the areas that you want to measure via the **"Point Selection Tool"**.
- Use the Zoom wheel in the smart panel below the computer monitor to zoom in (smaller square).
- Or out (larger square).


⚠ Be aware that the image pixel size and thus resolution changes accordingly.

Step 5 — Define your optimal imaging settings




- You can adjust your image settings through the **"Scanmode"** panel on the left and using the **"Live"**-button (gives you a preview with the actual imaging parameters).
- ☑ Proper xy sampling (pixel size) is crucial for acquiring optimal images.
- To adjust for the correct **pixel size** you can either use the online calculator such as the [SVI Nyquist Calculator](#) or the auto button for an estimate.
- Make sure laser intensity, line/frame accumulation/averaging are appropriate for the new format. Once you are happy you can start your image acquisition by pressing **"Start"**.

Step 6 — Optional - Define Z-stacks for different areas




- Define your **Z-stack** settings.
 - Use the smart wheel on the panel to go to the top of your sample and press **"Begin"**.
 - Go to the bottom and press **"End"**.
- If you want to define the Z-stacks individually for each area, make sure the **Same Stack Size for all Regions** is **"Off"**.
 - Go to the regions by double-clicking and repeat the steps above.
 - After defining beginning and end, click **"Redefine Stack"**.
 - A symbol appears indicating the redefined stack.
- Now you can start your image acquisition by pressing **"Start"**.

Step 7 — Optional - Tile Scans: Focus Map



- ❗ For more information on tile scans, please refer to the guide ["Leica LAS X Navigator"](#).
- If your sample is not 100% flat, it can happen that for tile scans not all the parts are in focus.
- Click on **"F+"** and add a series of focus points.
- ❗ Try to spread the focus points. The more you set, the better your sample will be in focus. However, do not set more than one focus point per tile.

Step 8 — Optional - Tile Scans: Focus Map



- Open **"Focus Map"** and double-click on the first focus point.
- Go to **"Fast Live"**.
- **Focus** in the respective field of view with the Z-wheel on the pannel.
- Click **"Set Z"**, which will save a Z-position for the focus point.
- Repeat these steps for all the focus points.
- When all focus points are set, close the window.
- Start your image acquisition by pressing **"Start"**.
 - All Z-positions of your tiles are now interpolated based on your focus points, thereby rendering all areas in focus.