

Leica Stellaris 5 inverse (USZ) - 1: Start-up

How to start up and mount a sample on the Leica Stellaris 5 inverse confocal laser scanning microscope located at USZ, Schmelzbergstrasse 12 , 8091 Zürich.

Written By: Joana Raquel Delgado Martins



INTRODUCTION



In this guide of the Center for Microscopy and Image Analysis we show how to start up and mount a sample on the NEW Leica Stellaris 5 inverse confocal laser scanning microscope located at the USZ, Schmelzbergstrasse 12 , CH-8091 Zürich.

It will introduce users in the new Stellaris 5 confocal laser scanning microscope. The most relevant new features are:

- NEW Spectrally flexible white light laser (WLL) with an extended spectral output on the red and near-infrared spectrum (485 nm up to 685 nm).
- NEW 3 Power HyD S detectors.
- NEW TauSense technology, which gives access to lifetime-based information, delivering additional insights to your experiments.


Please find more detailed information about the system [here](#).

Step 1 — Switching ON hardware




- Switch ON the **"Power"**, **"Laser"**, and turn the **"Emission"** key to "On-1" (control unit next to the microscope).
- The LED (for fluorescence observation via the oculars) switches on automatically. No action needed.
- Switch ON the PC.

Step 2 — Sign-In




- Sign-in with your **ZMB core** credentials.

Step 3 — Starting "LAS X" software





- **Start** the "LAS X" software.
 - Make sure "**machine.xlhw**" is selected as "Configuration", and "**DM18**" as "Microscope".
 - Click "**OK**".
 - Click "**Yes**" in order to **initialize the x/y stage**. *Please make sure nothing is placed currently on the stage.*
- ① *x/y stage initialization is necessary to be able to use the "LAS X Navigator".*

Step 4 — Choosing an objective



- **Lower the objective** by keeping the "LOWER Z" button pressed on the right side of the microscope.
- ⚠ This is a mandatory step as it avoids possible collision of the stage and objective during exchange of inserts and/or samples.
- In the **"Acquire" tab** you can now toggle between objectives (drop-down menu).
- Select the **5x dry objective** for easy sample navigation and focusing.
- ① *In order to facilitate the focusing it is recommended to start with a low magnification dry objective.*


Step 5 — Check stage insert



⚠ Make sure that the stage insert is correctly inserted and flat.



- Here correct and flat.
- Here not inserted correctly (stage not flat and shaky).

Step 6 — Mount and position your sample




- Push the condenser arm of the microscope back.
- Insert your sample with the **coverslip facing down** .
 - ❗ *Adjust the variable clamping range and moveable brackets to properly fix your sample.*
- Move your sample under the objective with the help of the external controller "**Smart Move**".
 - Movement in y-direction.
 - Movement in x-direction.
 - Toggle between coarse movement "XY Fast" and slow movement "XY Precise".
- ☑ Bring back condenser arm to its straight position.

Step 7 — Focus your sample




- On the touch screen at the microscope stand choose the **light path tab**.
- Click **"FLUO"** and choose an appropriate "FLUO-Filtercube" : e.g. "DAPI".
- Open the **"IL -Shutter"** (if activated the dot is yellow).
- Look through the oculars and **focus your sample** by using:
 - The z-wheel on the **external controller ("Smart Move")**.
 - ① Moving objectives upwards (towards sample) turn z-wheels clockwise/away from you. Moving objectives downwards (away from sample) turn z-wheels counter-clockwise/towards you.
 - Toggle between "Z FINE" and "Z COARSE" directly on the "Smart Move".
- If you cannot see any signal, make sure that the light path is directed to the eye pieces.

Step 8 — Optional - Save your focus position




- ① *The storage of the focal plane is helpful in order to find the focus back if the sample or objective will be changed.*
- To save your current focus position select the **"xyz tab"** and the **"Focusdrive Z"** on the touchscreen of the microscope.
- Click the **"Upper Focus Limit"** button.
- Press **"Set"**.
 - If done successfully you will see an upper marker line appearing.
- Press the **"Lower Limit"** button in order to move down (*for safe change of the objective or the sample*).

Step 9 — Switching to a higher magnification



- Lower the objectives, remove your sample and toggle within the software to the objective of choice.
- ① Depending on the objective different **immersion media** will be used. Apply directly on the sample.
 - Oil objectives: "Type-F" immersion liquid.
 - "Glycerin" objectives: "Type-G" immersion liquid.
 - "Water" objectives: ddH₂O (always use fresh).
- **Mount** your sample again and press the **"Upper Focus Limit"** button.
- ☑ Remember, you can move (back and forth) the condenser arm for ease of access.
- **Focus** your sample as described previously.

Step 10 — Additional information - Immersion objectives



- ⚠ For optimal imaging performance on some objectives the **correction collar** has to be adjusted.
 - **20x IMM** (multi-immersion - Oil, Glycerin or Water) needs to be set by moving the black dot to the corresponding **immersion media** ("OIL", "GLYC" or "0.17-W" (with cover glass) or "W-0" (without cover glass)).
 - Make sure that the cap of the **spring-loaded front lens** is released (working position). (*Please note: Displayed objective just an example, not exactly same as installed on system*).
- ⚠ Please, **DO NOT** remove the objectives for adjustment. They can be also accessed on the system.

Step 11 — Next steps



- You are now ready to acquire your images. Please check the corresponding follow up guides.
- [Leica Stellaris 5 - 2: Acquire an image in xy and z](#)